

EASY 2 CLEAN™

Silicone Full Face Mask



PneumoCare's Easy2Clean Mask is a Silicone Full Face Mask with easy to clean technology for use with ICU Ventilators or BiPAP/CPAP machines.

The unique mask design comes with an Auto Adjustable Frame and a Fused Dual Wall Soft Silicone Cushioning which not only provides exceptional comfort to the patient but also delivers in terms of quality and performance.

Mask Categories



Non-Vented
Easy2Clean Non-Vented mask is designed to be used with double limb circuit ICU Ventilators

Vented
Easy2Clean Vented mask is designed to be used with single limb circuit BiPAP/CPAP machines

