


Easy2Neb™

NIV Nebulizer Full Face Mask



PneumoCare's Easy2Neb Mask is world's first NIV Full Face Mask with universal Nebulization port.

- Easy2Neb Mask offers Aerosol Delivery during Non-Invasive Ventilation
- The Nebulizer port is between the expiratory port and the mask leading to higher deposition rates

The unique mask design comes with an Auto Adjustable Frame and a Fused Dual Wall Soft Silicone Cushioning which not only provides exceptional comfort to the patient but also delivers in terms of quality and performance.

Mask Categories



Vented

Easy2Neb Vented mask is designed to be used with single limb circuit BiPAP/CPAP machines

Non-Vented
Easy2Neb Non-Vented mask is designed to be used with double limb circuit ICU Ventilators

