

Tracheotomy Tube

Cuff (Normal Type)

Cat.No	Size
TYN0025	2.5
TYN0030	3.0
TYN0035	3.5
TYN0040	4.0
TYN0045	4.5
TYN0050	5.0
TYN0055	5.5
TYN0060	6.0
TYN0065	6.5

Cat.No	Size
TYN0070	7.0
TYN0075	7.5
TYN0080	8.0
TYN0085	8.5
TYN0090	9.0
TYN0095	9.5
TYN0100	10
TYN0110	11.

21



Features

- ◆ To establish artificial respiration channels.
- ◆ Movable connecting ring embedded for good compliance.
- ◆ The patient-end of the obturator is bullet-shaped, convenient for inserting into the incision.
- ◆ Radiopaque line (X-ray visualization).

Uncuff (Normal Type)

Cat.No	Size
TYN1025	2.5
TYN1030	3.0
TYN1035	3.5
TYN1040	4.0
TYN1045	4.5
TYN1050	5.0
TYN1055	5.5
TYN1060	6.0
TYN1065	6.5

Cat.No	Size
TYN1070	7.0
TYN1075	7.5
TYN1080	8.0
TYN1085	8.5
TYN1090	9.0
TYN1095	9.5
TYN1100	10.0
TYN1110	11.0



Features

- ◆ To establish artificial respiration channels.
- ◆ Movable connecting ring embedded for good compliance.
- ◆ The patient-end of the obturator is bullet-shaped, convenient for inserting into the incision.
- ◆ Radiopaque line (X-ray visualization).
- ◆ Uncuffed tube body.