

OXY SET TWIN BORE NASAL SET



14

OXY SET® GS-2016

- Oxy Set is suitable for easy application and efficient administration of oxygen concentration.
- The system of attachment provides maximum freedom to the patient and leaves the patient's mouth free for nutritional and communication purposes.
- Twin prong/ nasal tips are designed to ensure equal volume of oxygen to both the air passages, of the patient.
- Soft funnel shaped connector facilitates easy connection to oxygen source.
- Multichannel tube ensures continuous supply of oxygen, even if the tube is accidentally kinked.
- Manufactured from soft, light weight non toxic material.
- **AVAILABLE OPTIONS :**
- Sterile : Peelable Pack ; Box of 20/Master box of 200.
- Non Sterile : Packed in polybag ; Master box of 100.

SPECIFICATIONS

Prong Size	Tube Length
Adult	210 cm
Child	210 cm
Infant	210 cm
Neonate	210 cm