



## **NAREENA'S PORTABLE SUCTION MACHINE**

(Model No: NLS- PSM001)

Portable suction machine is a new design with an oil-free lubricated pump, developed after a long research and development. It is specially designed to suck thick liquid such as blood, phlegm or any other liquid waste. The product is widely accepted and used across hospitals, medical institutions and homes. The machine is crafted with high-grade plastic parts for better finish and durability of the product. Silent pump maintains low noise-level while maintaining proper pressure to suck liquid waste. Fitted with in plastic 1000 ml reservoir.

### Specifications:

Power supply: AC220V $\pm$ 10% 50HZ $\pm$ 2%

Max negative pressure	:	$\geq$ 0.075Mpa
Noise	:	$\geq$ 60dB(A)
Power	:	90VA
Pumping rate	:	$\geq$ 18L/min
Reservoir capacity	:	> 1000ml"

