

# GUEDEL AIRWAYS – Oropharyngeal Airway



## PRODUCT CODES

CODE	SIZE	LENGTH	COLOUR
SAHGUD000	Size 000	40mm	Pink
SAHGUD001	Size 00	50mm	Blue
SAHGUD002	Size 0	60mm	Black
SAHGUD003	Size 1	70mm	White
SAHGUD004	Size 2	80mm	Green
SAHGUD005	Size 3	90mm	Yellow
SAHGUD006	Size 4	100mm	Red
SAHGUD007	Size 5	110mm	L. Blue

- An oropharyngeal airway (also known as a Guedel airway) is an airway adjunct used to maintain a patient's airway.
- Suitable for maintaining an unobstructed oropharyngeal airway during general anaesthesia and in unconscious patients.
- Commonly used in pre-hospital emergency care and for short term airway management.

## FEATURES AND BENEFITS

- **Rounded atraumatic edges allows for easy insertion**
- **Smooth airway path for easy cleaning**
- **Bite block to prevent biting of tongue and airway occlusion**
- **Colour coded for easy identification of sizes – as per ISO5364 specifications**
- **Available in sizes: 000, 00, 0, 1, 2, 3, 4 & 5**
- **Individually packaged sterile product for single use**

## SPECIFICATIONS & MATERIALS

<b>Guedel Airways</b>	- <b>100% PVC</b>
<b>Sterility</b>	- <b>ETO</b>
	- <b>Latex Free</b>



STERILE EO



South African  
Hospital Supplies Pty. Ltd.

P.O. Box 26852 Hout Bay South Africa 7872 Tel: +27 21 671 0796 Fax: +27 21 671 1189 www.sahospitalssupplies.co.za

## BRANCH CONTACT DETAILS

Head Office – Cape Town: Tel: 021-671 0796

Port Elizabeth: 084 609 5936

KwaZulu Natal: 074 677 1122

Johannesburg: 011-450 1158