

Omron Blood Pressure Monitor HBP-1320



OM-HBP1320

Product Features

- Dual Measurement Mode
- Motion stop function
- Irregular pulse Indicator (disabled when used with auscultation mode)
- Helps identify changes in heart rate, rhythm, or pulse which may be caused by heart disease or other serious health issues.
- Last reading display
- Long Battery life
- Custom rechargeable battery makes it possible to provide 300 measurements with one charge.
- Power saving Auto Power Off
- 5 cuff sizes available for measurement
- the unit can be used on pediatrics of age 3 years and above.
- Cuff sizes available for use on different arm circumference: XL: 42-50cm, M: 22-32cm, S:17-22cm, SS:12-18cm