


Spiroflo® Filter

For Spirometry Use



Description

Disposable filters for pulmonary function test with electrostatically charged filter medium for high efficiency filtration of BEF is 99.99990% and VEF is 99.99975%. Can be incorporated with patient friendly mouthpiece with negligible resistance for accurate measurements. Available in combo packs with mouthpiece and nose clip.

Specification

- Housing: White polypropylene
- Electrostatic charged synthetic fibre polypropylene non-woven media
- BEF is 99.99990% & VEF is 99.99975%
- Weight: 27gms
- Height: 8.5cm
- Resistance < 2.3 mm H₂O @ 32 LPM
- Dead volume 75ml
- Latex Free

Order Information

Cat No	Description	MOQ
5501	Spiroflo Filter (Oval)	100 nos
5511	Spiroflo Filter (Round)	100 nos



JK Medical Systems Private Limited

(AN ISO 13485:2016 & WHO-GMP CERTIFIED COMPANY)

1068, Munusamy Salai, West KK Nagar, Chennai-600 078, India | Ph:044-23664761/ 62 | sales@jkmedical.co.in
www.jkmedical.co.in