

25E IND ALGAN S0002  
84MM 9.00KG 10L PW150  
EN ISO 7866:2012



**MEDICAL/INDUSTRIAL/  
SPECIALTY GASES CYLINDERS**



**SCUBA CYLINDERS**



**FIRE EXTINGUISHER  
CYLINDERS**



**BEVERAGES CYLINDERS**

**Aluminium High  
Pressure Gas Cylinders ...**

# Oxy<sup>kit</sup><sup>®</sup>


Home Care



Easy To Operate




**IN EMERGENCY**



**RESPIRATORY PROBLEM**



**HIGH ALTITUDE & TREKKING**



**EASY TRAVEL**

OxyKit comes with Al-Can make Light Weight Aluminium Oxygen Cylinder with Valve, Regulator (0-15L Flow rate adjustable), Mask and Carry bag. The Cylinder and Valve are Approved by Government of India - PESO (Petroleum and Explosive Department). OxyKit Aluminium Cylinder is approx 30% Light Weight compared to Steel Cylinder. Also it is Portable and Easy to use. The Cylinder is Refillable which can be Refilled anywhere in India from any Gas Manufacturers or Gas Dealers.

#### OXYKIT IS IDEAL FOR :

1. Patients suffering from Respiratory Disorder (under Doctor's guidance/prescription).
2. Senior Citizens for Well-Being.
3. Stress Relieving.
4. Yatra, Mountaineers, Trekker for Ease of Breathing at High Altitude.
5. Athletes and Sportsmen.
6. Refresh from daily Pollution.



**1.8 ltr**



**2.2 ltr**



**3.1 ltr**



**4.5 ltr**



**9.1 ltr**



**10 ltr**



**12 ltr**



**LIGHT WEIGHT**



# Oxy<sup>kit</sup><sup>®</sup>

## Hospital & Ambulance



Corrosion Free



N2O & Co2 CYLINDER FOR ANESTHESIA MACHINE



O2 CYLINDER FOR BABY INCUBATOR



Co2 CYLINDER FOR LAPAROSCOPY MACHINE



HELIUM CYLINDER FOR IABP MACHINE

Oxykit's Lightweight, Portable & Non Magnetic Aluminium Cylinders are widely used for Medical Gases in ICU, NICU, Operation Theatre, Emergency Ward, Incubator, Anesthesia Machine, Stretcher and Ambulances



100L OxyKit can be used for Symptoms of Headache, Helping in wound Healing, Senior Citizen needing Supplement oxygen, A useful Kit for Mountaineers, Trekkers, Sports Person and at the High Altitude to Breathe easy.



Oxynox, Helps To Give Instant Energy To Boost Your Power Anytime... Anywhere.  
Gas Content: 9ltr.@25°C (Appox.)  
Weight: 140Gms with gas ( Approx.)  
Net Volume: 498ml (Approx.)



Mixture of Nitrous Oxide & Oxygen Gas Kit is also known as Pain Relief Gas has been used in Pain Relief, Sedation and Anesthesia, Since many years.



Oxygen (O<sub>2</sub>)



ENTONOX  
(O<sub>2</sub> + N<sub>2</sub>O)



CARBON DIOXIDE  
(CO<sub>2</sub>)



NITROUS OXIDE  
(N<sub>2</sub>O)



AIR



PORTABLE

Approved By Petroleum And Explosive Safety Organisation (PESO), Govt Of India.

ISO 7866 With "PI" (π) Mark (TPED)/ CE Mark (PED) Approved By Apragaz.



PESO APPROVED



CE MARK (PED)  
APPROVED



PIE MARK (TPED)  
APPROVED



DIMENSION CHECK



LEAK TEST



HYDRAULIC TEST



INTERNAL SURFACE  
CLEANLINESS CHECK

## Medical/Industrial Cylinders Working Pressure 150 Bar

Item No.	Water Capacity	Outside Dia	Gas Capacity		Length And Weight	
	Liters	mm	C.ft	Liters	Inch	Kgs
ACE-0.495L-65	0.495	65	3	75	9	0.5
ACE-1L-90	1	90	5	150	9	1.2
ACE-1.5L-90	1.5	90	8	225	14	1.8
ACE-1.8L-111	1.8	111	10	270	11	1.9
ACE-2.2L-111	2.2	111	12	330	13	2.2
ACE-3.1L-111	3.1	111	16	465	18	3
ACE-4.5L-111	4.5	111	24	675	25	4.1
ACE-5L-140	5	140	27	750	19	4.9
ACE-6.8L-140	6.8	140	36	1020	25	6.4
ACE-7.8L-140	7.8	140	41	1170	28	7.1
ACE-9.1L-140	9.1	140	48	1365	32	7.6
ACE-10L-152	10	152	53	1500	32	10.1
ACE-10L-176	10	176	53	1500	23	9.3
ACE-12L-176	12	176	64	1800	27	10.8
ACE-15L-176	15	176	79	2250	33	13.1
ACE-20L-203	20	203	106	3000	33	16.9
ACE-21L-203	21	203	111	3150	35	17.6

# LIFE SAVER!

When every second counts



NON MAGNETIC



# Oxy<sup>kit</sup><sup>®</sup>

Accessories

## Valves

**Bull Nose Valve**



Wheel Type



**Pin Index Valve**



Toggle Type



## Regulators

**Bull Nose Regulator**



**Pin Index Regulator**



## Valve Guard

**Metal/  
Aluminium  
Valve Guard**



**Plastic  
Valve Guard**



## Trolley & Bag

**Trolley**



**Bag**





## Aluminum High Pressure Gas Cylinders ...

### **AL-CAN EXPORTS PVT. LTD.**

Sheetal Indl. Estate, Mira-Bhayander Road,  
Bhayander East, Dist.- Thane. Pin - 401 105,  
Tel : +91 22 2814 2478  
Email: [cylinder@alcanexports.com](mailto:cylinder@alcanexports.com)  
Web Site: [www.alcancylinder.com](http://www.alcancylinder.com)

---

### **AL-CAN INTERNATIONAL LLC**

P.O. BOX: 84685, Al Qusais Ind. Area-2, Dubai, UAE  
Tel: +971 4 2671101 | Fax: +971 4 2582490  
Email id: [info@alcaninternational.com](mailto:info@alcaninternational.com) |  
Web Site: [www.alcaninternational.com](http://www.alcaninternational.com) | [www.alcancylinder.com](http://www.alcancylinder.com)