



**MEDICAL/INDUSTRIAL/
SPECIALITY GASES CYLINDERS**



SCUBA CYLINDERS



**FIRE EXTINGUISHER
CYLINDERS**



BEVERAGES CYLINDERS

**Aluminium High
Pressure Gas Cylinders ...**

Oxykit[®]

Home Care




Easy To Operate



IN EMERGENCY



HIGH ALTITUDE & TREKKING



RESPIRATORY PROBLEM



EASY TRAVEL

OxyKit comes with Al-Can make Light Weight Aluminium Oxygen Cylinder with Valve, Regulator (0-15L Flow rate adjustable), Mask and Carry bag. The Cylinder and Valve are Approved by Government of India - PESO (Petroleum and Explosive Department). OxyKit Aluminium Cylinder is approx 30% Light Weight compared to Steel Cylinder. Also it is Portable and Easy to use. The Cylinder is Refillable which can be Refilled anywhere in India from any Gas Manufacturers or Gas Dealers.

OXYKIT IS IDEAL FOR :

1. Patients suffering from Respiratory Disorder (under Doctor's guidance/prescription).
2. Senior Citizens for Well-Being.
3. Stress Relieving.
4. Yatra, Mountaineers, Trekker for Ease of Breathing at High Altitude.
5. Athletes and Sportsmen.
6. Refresh from daily Pollution.



1.8 ltr

2.2 ltr

3.1 ltr

4.5 ltr

9.1 ltr

10 ltr

12 ltr

 LIGHT WEIGHT

Oxykit[®]

Hospital
& Ambulance



Corrosion Free



N2O & Co2 CYLINDER FOR ANESTHESIA MACHINE



O2 CYLINDER FOR BABY INCUBATOR



Co2 CYLINDER FOR LAPAROSCOPY MACHINE



HELIUM CYLINDER FOR IABP MACHINE

Oxykit's Lightweight, Portable & Non Magnetic Aluminium Cylinders are widely used for Medical Gases in ICU, NICU, Operation Theatre, Emergency Ward, Incubator, Anesthesia Machine, Stretcher and Ambulances



100L OxyKit can be used for Symptoms of Headache, Helping in wound Healing, Senior Citizen needing Supplement oxygen, A useful Kit for Mountaineers, Trekkers, Sports Person and at the High Altitude to Breathe easy.



Oxynox, Helps To Give Instant Energy To Boost Your Power Anytime... Anywhere.
Gas Content: 9ltr.@25°C (Approx.)
Weight: 140Gms with gas (Approx.)
Net Volume: 498ml (Approx.)



Mixture of Nitrous Oxide & Oxygen Gas Kit is also known as Pain Relief Gas has been used in Pain Relief, Sedation and Anesthesia, Since many years.



Oxygen (O₂)



ENTONOX
(O₂ + N₂O)



CARBON DIOXIDE
(CO₂)



NITROUS OXIDE
(N₂O)



AIR



PORTABLE

Approved By Petroleum And Explosive Safety Organisation (PESO), Govt Of India.
ISO 7866 With "PI" (π) Mark (TPED)/ CE Mark (PED) Approved By Apragaz.



PESO APPROVED



CE MARK (PED)
APPROVED



PIE MARK (TPED)
APPROVED



DIMENSION CHECK



LEAK TEST



HYDRAULIC TEST



INTERNAL SURFACE
CLEANLINESS CHECK

Medical/Industrial Cylinders Working Pressure 150 Bar

Item No.	Water Capacity	Outside Dia	Gas Capacity		Length And Weight	
	Liters		C.ft	Liters	Inch	Kgs
ACE-0.495L-65	0.495	65	3	75	9	0.5
ACE-1L-90	1	90	5	150	9	1.2
ACE-1.5L-90	1.5	90	8	225	14	1.8
ACE-1.8L-111	1.8	111	10	270	11	1.9
ACE-2.2L-111	2.2	111	12	330	13	2.2
ACE-3.1L-111	3.1	111	16	465	18	3
ACE-4.5L-111	4.5	111	24	675	25	4.1
ACE-5L-140	5	140	27	750	19	4.9
ACE-6.8L-140	6.8	140	36	1020	25	6.4
ACE-7.8L-140	7.8	140	41	1170	28	7.1
ACE-9.1L-140	9.1	140	48	1365	32	7.6
ACE-10L-152	10	152	53	1500	32	10.1
ACE-10L-176	10	176	53	1500	23	9.3
ACE-12L-176	12	176	64	1800	27	10.8
ACE-15L-176	15	176	79	2250	33	13.1
ACE-20L-203	20	203	106	3000	33	16.9
ACE-21L-203	21	203	111	3150	35	17.6

LIFE SAVER!

When every second counts



NON MAGNETIC

Oxy *kit*[®] Accessories

Valves

Bull Nose Valve



Wheel Type



Pin Index Valve



Toggle Type




Regulators

Bull Nose Regulator



Pin Index Regulator



Valve Guard

Metal/Aluminium Valve Guard



Plastic Valve Guard



Trolley & Bag

Trolley



Bag





Aluminum High Pressure Gas Cylinders ...

AL-CAN EXPORTS PVT. LTD.

Sheetal Indl. Estate, Mira-Bhayander Road,
Bhayander East, Dist.- Thane. Pin - 401 105,
Tel : +91 22 2814 2478

Email: cylinder@alcanexports.com

Web Site: www.alcancylinder.com

AL-CAN INTERNATIONAL LLC

P.O. BOX: 84685, Al Qusais Ind. Area-2, Dubai, UAE

Tel: +971 4 2671101 | Fax: +971 4 2582490

Email id: info@alcaninternational.com |

Web Site: www.alcaninternational.com | www.alcancylinder.com